

Mouse Squamous Cell Carcinoma Antigen (SCCA) ELISA Kit



Catalog No: #EK20788

Orders: order@signalwayantibody.com

Package Size: #EK20788-1 48T #EK20788-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Squamous Cell Carcinoma Antigen (SCCA) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:62.5-4000pg/mL

Sensitivity:23.6pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only