

# Porcine Beta-2-Microglobulin (b2M) ELISA Kit

Catalog No: #EK20789



Package Size: #EK20789-1 48T #EK20789-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Porcine Beta-2-Microglobulin (b2M) ELISA Kit |
| Brief Description  | ELISA Kit                                    |
| Applications       | ELISA  |
| Species Reactivity | Porcine                                      |

## Application Details

Detect Range:3.12-200ng/mL

Sensitivity:1.27ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only