

Recombinant Human 5-hydroxytryptamine receptor 1A (HTR1A)



Catalog No: #AP89688

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89688-1 20ug #AP89688-2 100ug #AP89688-3 1mg

Description

Product Name	Recombinant Human 5-hydroxytryptamine receptor 1A (HTR1A)
Brief Description	Recombinant Protein
Other Names	HTR1A; ADRB2RL1; ADRBRL1; 5-hydroxytryptamine receptor 1A; 5-HT-1A; 5-HT1A; G-21; Serotonin receptor 1A
Uniprot	P08908
Target Species	Homo sapiens (Human)
Tag Info	Tag type will be determined during the manufacturing process.
Target Sequence	Protein Length:Partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Application Details

>85% (SDS-PAGE)

Note: This product is for in vitro research use only