

## Anti-CD86 Antibody

Catalog No: #62034



Package Size: #62034-1 50ul #62034-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Anti-CD86 Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1 kappa
Purification	The antibody was purified by affinity chromatography.
Applications	IF, FC
Species Reactivity	H
Immunogen Description	Native purified human CD86.
Target Name	CD86
Other Names	CD28LG2; T-lymphocyte activation antigen CD86; Activation B7-2 antigen; B70; BU63; CTLA-4 counter-receptor B7.2; FUN-1; CD antigen CD86
Uniprot	P42081
GeneID	942
Formulation	Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note: This product is for in vitro research use only