

CD206 mouse mAb

Catalog No: #49238



Package Size: #49238-1 50ul #49238-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD206 mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1, Kappa
Purification	Protein G
Applications	WB;IF;ELISA
Species Reactivity	Human;Mouse;Rat;
Specificity	This antibody detects endogenous levels of CD206.
Immunogen Description	Synthesized peptide derived from human protein.
Target Name	MRC1
Other Names	MRC1;CLEC13D;CLEC13DL;MRC1L1;Macrophage mannose receptor 1;MMR;C-type lectin domain family 13 member D;C-type lectin domain family 13 member D-like;Macrophage mannose receptor 1-like protein 1;CD206
Uniprot	P22897
GeneID	4360
Calculated MW	190kD
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000. IF 1:100-500. ELISA 1:1000-5000

Note: This product is for in vitro research use only