

## PDX1 Rabbit pAb

Catalog No: #62004



Package Size: #62004-1 50ul #62004-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PDX1 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	IHC-P, IHC-F, Flow-Cyt, IF
Species Reactivity	Human, Mouse, Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human PDX1: 101-200/284
Other Names	Glucose sensitive factor; Glucose-sensitive factor; GSF; IDX 1; IDX-1; IDX1; Insulin promoter factor 1; insulin promoter factor 1 homeodomain transcription factor; insulin upstream factor 1; IPF 1; IPF-1; IPF1; Islet/duodenum homeobox 1; Islet/duodenum homeobox-1; IUF 1; IUF-1; IUF1; MODY4; Pancreas/duodenum homeobox 1; Pancreas/duodenum homeobox protein 1; pancreatic and duodenal homeobox P; PDX 1; PDX-1; PDX1_HUMAN; Somatostatin transactivating factor 1; Somatostatin-transactivating factor 1; STF 1; STF-1; STF1.
Uniprot	P52946
GeneID	3651
Calculated MW	30 kDa
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P 1:100-500 IHC-F 1:100-500 Flow-Cyt 1&amp;mu;g/Test IF 1:100-500

Note: This product is for in vitro research use only