

LIAS rabbit pAb

Catalog No: #62005



Package Size: #62005-1 50ul #62005-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LIAS rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;WB
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of LIAS at Human, Mouse,Rat
Immunogen Description	Synthesized peptide derived from human N-terminal LIAS
Target Name	LIAS LAS HUSSY-01
Other Names	Lipoyl synthase, mitochondrial (EC 2.8.1.8) (Lipoate synthase) (LS) (Lip-syn) (Lipoic acid synthase)
Uniprot	O43766
GeneID	11019
Calculated MW	41kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000 IHC 1:50-200

Note: This product is for in vitro research use only