

ULK1 Polyclonal Antibody

Catalog No: #62009

Package Size: #62009-1 50ul #62009-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ULK1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;IF;WB;ELISA
Species Reactivity	Human;Mouse
Specificity	ULK1 Polyclonal Antibody detects endogenous levels of ULK1 protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human ULK1.
Target Name	ULK1
Other Names	ULK1;KIAA0722;Serine/threonine-protein kinase ULK1;Autophagy-related protein 1 homolog;ATG1;hATG1;Unc-51-like kinase 1
Uniprot	O75385
GeneID	8408
Calculated MW	113kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000,IHC 1:100 - 1:300. ELISA: 1:40000.. IF 1:50-200

Images



Immunohistochemical analysis of paraffin-embedded human tonsil.

Note: This product is for in vitro research use only