

Human Cathepsin L (CTSL) ELISA Kit

Catalog No: #EK20795



Package Size: #EK20795-1 48T #EK20795-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Cathepsin L (CTSL) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	CTSL1; CATL; CTS-L; MEP; Cathepsin L1; Major excreted protein

Application Details

Detect Range:62.5-4000pg/mL

Sensitivity:21.9pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only