

## F-Actin Rabbit pAb

Catalog No: #62036



Package Size: #62036-1 50ul #62036-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	F-Actin Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC- Flow-Cyt IF
Species Reactivity	Human, Mouse, Rat
Immunogen Description	Synthetic MAP peptide derived from human beta-Actin: 1-50/374
Other Names	f actin; Filamentous actin; F-actin capping protein alpha subunit; CapZ alpha-1; CAZA1_HUMAN; CAPZA1.
Uniprot	P60709
GeneID	60
Calculated MW	42 kDa
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:1000-5000, IHC-P=1:100-500, IHC-F=1:100-500, Flow-Cyt=3&amp;mu;g /test, IF=1:100-500

Note: This product is for in vitro research use only