

Human SPARC-related modular calcium-binding protein 1 (SMOC1) ELISA Kit



Catalog No: #EK20751

Orders: order@signalwayantibody.com

Package Size: #EK20751-1 48T #EK20751-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human SPARC-related modular calcium-binding protein 1 (SMOC1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.054ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only