

NECD rabbit pAb

Catalog No: #62040



Package Size: #62040-1 50ul #62040-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NECD rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of human NECD
Immunogen Description	Synthesized peptide derived from human NECD
Target Name	NDN
Uniprot	Q99608
GeneID	4692
Calculated MW	35kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only