

# ITGB5 Mouse Monoclonal Antibody

Catalog No: #62041

Package Size: #62041-1 50ul #62041-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	ITGB5 Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IF,FCM
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human ITGB5 (AA: 50-319) expressed in E. Coli.
Uniprot	P18084
GeneID	3693
Calculated MW	88kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details


IF 1/50 - 1/200

FCM 1/200 - 1/400

## Images



Immunofluorescence analysis of HeLa cells using ITGB5 mouse mAb (green).



Flow cytometric analysis of Jurkat cells using ITGB5 mouse mAb (green) and negative control (red).

---

Note: This product is for in vitro research use only