

Necdin Rabbit pAb

Catalog No: #62044



Package Size: #62044-1 50ul #62044-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Necdin Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F IF
Species Reactivity	Human, Mouse, Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human Necdin: 101-200/321
Other Names	HsT16328; NDN; NECD_HUMAN; Necdin; Necdin homolog.
Uniprot	Q99608
GeneID	4692
Calculated MW	36 kDa
Concentration	1 mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:100-500 IHC-F=1:100-500 IF=1:100-500

Note: This product is for in vitro research use only