

Mouse Nociceptin receptor α OPRL1 α ELISA Kit

Catalog No: #EK20801



Package Size: #EK20801-1 48T #EK20801-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Nociceptin receptor α OPRL1 α ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:78-5000pg/mL
Sensitivity:37.5pg/mL
Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids
Sample Volume: 50 μ L
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only