

BAI1 Polyclonal Antibody

Catalog No: #62045

Package Size: #62045-1 50ul #62045-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	BAI1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB,IHC,IF
Species Reactivity	Human,Rat,Mouse
Specificity	The antibody detects endogenous BAI1 protein
Immunogen Description	Synthetic Peptide of BAI1
Target Name	BAI1
Other Names	Brain-specific angiogenesis inhibitor 1
Uniprot	O14514
GeneID	575
Calculated MW	130-170kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000. IHC 1:100 - 1:300. IF 1:50-200

Images



Western blot analysis of 1) 293T, 2) Mouse Brain Tissue, 3) Rat Brain Tissue with BAI1 Rabbit pAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using BAI1 Rabbit pAb diluted at 1:200.

Note: This product is for in vitro research use only