

# Human RNA-binding protein with serine-rich domain 1 ELISA Kit



Catalog No: #EK20802

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20802-1 48T #EK20802-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human RNA-binding protein with serine-rich domain 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:13.5pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only