

Human 39S ribosomal protein L21, mitochondrial (MRPL21) ELISA Kit



Catalog No: #EK20809

Orders: order@signalwayantibody.com

Package Size: #EK20809-1 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human 39S ribosomal protein L21, mitochondrial (MRPL21) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	L21mt; MGC62013; MRP-L21; 39S ribosomal protein L21; mitochondrial
Uniprot	Q2TBS2

Application Details

Detect Range:Request Information
Sensitivity:Request Information
Sample Type:Serum, Plasma, Other biological fluids
Sample Volume: 1-200 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only