

# Human 39S ribosomal protein L23, mitochondrial (MRPL23) ELISA Kit



Catalog No: #EK20813

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20813-1 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human 39S ribosomal protein L23, mitochondrial (MRPL23) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	FLJ45387; L23MRP; RPL23; RPL23L; L23 mitochondrial-related protein OTTHUMP00000069610 ribosomal protein L23-like ribosomal protein related to L23 (mitochondrial)
Uniprot	Q16540

## Application Details

Detect Range:Request Information
Sensitivity:Request Information
Sample Type:Serum, Plasma, Other biological fluids
Sample Volume: 1-200 &mu;L
Assay Time:1-4.5h
Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only