

Human 28S ribosomal protein S29, mitochondrial (DAP3) ELISA Kit



Catalog No: #EK20821

Orders: order@signalwayantibody.com

Package Size: #EK20821-1 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human 28S ribosomal protein S29, mitochondrial (DAP3) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	DAP-3; DKFZp686G12159; MGC126058; MGC126059; MRP-S29; MRPS29; bMRP-10; death-associated protein 3 mitochondrial 28S ribosomal protein S29
Uniprot	P51398

Application Details

Detect Range:Request Information
Sensitivity:Request Information
Sample Type:Serum, Plasma, Other biological fluids
Sample Volume: 1-200 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only