

CTR1 rabbit pAb

Catalog No: #62051

Package Size: #62051-1 50ul #62051-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CTR1 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	This antibody detects endogenous levels of CTR1 at Human/Mouse/Rat
Immunogen Description	Synthesized peptide derived from human CTR1
Target Name	SLC7A1 ATRC1 ERR REC1L
Uniprot	P30825
GeneID	6541
Calculated MW	69kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1 μ g 500-2000

Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only