

Annexin V Rabbit pAb, APC conjugated

Catalog No: #C90188APC



Package Size: #C90188APC-1 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Annexin V Rabbit pAb, APC conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Human Mouse Rat
Immunogen Description	recombind human Annexin V
Conjugates	APC
Other Names	Anchorin CI; Annexin 5; Annexin A5; Annexin V; Annexin5; AnnexinA5; AnnexinV; ANX 5; ANX A5; ANX5; ANXA5; Calphobindin I; CBP I; CBP-I; Endonexin II; ENX 2; ENX2; Lipocortin V; PAP I; Placental anticoagulant protein I; PAP-I; Placental anticoagulant protein 4; PP 4; PP4; Thromboplastin inhibitor; VAC alpha; Vascular anticoagulant alpha; Vascular anticoagulant-alpha; ANXA5_HUMAN; Anchorin CII; Annexin-5; VAC-alpha.
Uniprot	P08758
GeneID	308
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1µg /test, IF=1:100-500

Note: This product is for in vitro research use only