

# Rat Plasmin-antiplasmin complex (PAP) ELISA Kit

Catalog No: #EK20805



Package Size: #EK20805-1 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Rat Plasmin-antiplasmin complex (PAP) ELISA Kit |
| Brief Description  | ELISA Kit                                       |
| Applications       | ELISA   |
| Species Reactivity | Rat   |

## Application Details

Detect Range:3.12-200ng/mL

Sensitivity:1.56ng/mL

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only