

Rat Plasmin-antiplasmin complex (PAP) ELISA Kit

Catalog No: #EK20805



Package Size: #EK20805-1 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Plasmin-antiplasmin complex (PAP) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:3.12-200ng/mL

Sensitivity:1.56ng/mL

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only