

# Human Intracellular hyaluronan-binding protein 4(HABP4) ELISA Kit



Catalog No: #EK20816

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20816-1 48T #EK20816-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Human Intracellular hyaluronan-binding protein 4(HABP4) ELISA Kit |
| Brief Description  | ELISA Kit   |
| Applications       | ELISA   |
| Species Reactivity | Human   |

## Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.116ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only