

Human Rap guanine nucleotide exchange factor 2(RAPGEF2) ELISA Kit



Catalog No: #EK20825

Orders: order@signalwayantibody.com

Package Size: #EK20825-1 48T #EK20825-2 96T

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Human Rap guanine nucleotide exchange factor 2(RAPGEF2) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Human |

Application Details

| |
|--|
| Detect Range:0.156-10ng/mL |
| Sensitivity:0.073ng/mL |
| Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids |
| Sample Volume: 50 μL |
| Assay Time:1-4.5h |
| Detection wavelength:450 nm |

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only