

PLVAP Rabbit pAb

Catalog No: #62052



Package Size: #62052-1 50ul #62052-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PLVAP Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC/IF
Species Reactivity	Human Mouse Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human PLVAP: 301-400/442
Other Names	FELS; Fenestrated endothelial linked structure protein; Fenestrated endothelial-linked structure protein; gp68; Hypothetical protein DKFZp667O074; Meca 32 antigen; MECA32 antigen; Panendothelial cell antigen; Plasmalemma vesicle associated protein; Plasmalemma vesicle protein 1; Plasmalemma vesicle-associated protein; Plvap; PLVAP_HUMAN; PV 1; PV-1; PV1; PV1 protein.
Uniprot	Q9BX97
GeneID	83483
Calculated MW	50kDa
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, ICC/IF=1:100-500

Note: This product is for in vitro research use only