

## Ribosomal Protein S21 Polyclonal Antibody

Catalog No: #62054

Package Size: #62054-1 50ul #62054-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Ribosomal Protein S21 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC IF ELISA
Species Reactivity	Human Mouse Rat
Specificity	Ribosomal Protein S21 Polyclonal Antibody detects endogenous levels of Ribosomal Protein S21 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human RPS21.
Target Name	RPS21
Other Names	RPS21 ; 40S ribosomal protein S21
Uniprot	P63220
GeneID	6227
Calculated MW	9kDa
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

IHC 1:100-1:300;ELISA 1:10000;IF 1:50-200

## Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using RPS21 Antibody. The picture on the right is blocked with the synthesized peptide.

---

Note: This product is for in vitro research use only