

## CUL9 rabbit pAb

Catalog No: #62055

Package Size: #62055-1 50ul #62055-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CUL9 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC WB
Species Reactivity	Human Mouse
Specificity	This antibody detects endogenous levels of CUL9 at Human, Mouse
Immunogen Description	Synthesized peptide derived from human C-terminal CUL9
Target Name	CUL9 H7AP1 KIAA0708 PARC
Other Names	Cullin-9 ; CUL-9 ; UbcH7-associated protein 1 ;
Uniprot	Q8IWT3
GeneID	23113
Calculated MW	277kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:500-2000;IHC 1:50-200

Note: This product is for in vitro research use only