

## PANK3 rabbit pAb

Catalog No: #62057



Package Size: #62057-1 50ul #62057-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PANK3 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human Mouse
Specificity	This antibody detects endogenous levels of human PANK3
Immunogen Description	Synthesized peptide derived from human PANK3
Target Name	PANK3
Uniprot	Q9H999
GeneID	79646
Calculated MW	41kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only