

TNNT1 rabbit polyclonal antibody

Catalog No: #62059



Package Size: #62059-1 50ul #62059-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	TNNT1 rabbit polyclonal antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC WB
Species Reactivity	Human Mouse
Immunogen Description	Fusion protein of human TNNT1
Other Names	troponin T1, slow skeletal type; ANM; TNT; NEM5; STNT; TNTS
Uniprot	P13805
GeneID	7138
Concentration	0.7 mg/ml
Formulation	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at -20Λ C

Application Details

WB 1:500-2000;IHC 1:25-50

Note: This product is for in vitro research use only