

Rat WNT1 Inducible Signaling Pathway Protein 1 (WISP1) ELISA Kit

Catalog No: #EK20820

Orders: order@signalwayantibody.com

Package Size: #EK20820-1 48T #EK20820-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Rat WNT1 Inducible Signaling Pathway Protein 1 (WISP1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:6.3pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only