

Mouse Forkhead box protein P1(FOXP1) ELISA Kit

Catalog No: #EK20830



Package Size: #EK20830-1 48T #EK20830-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Forkhead box protein P1(FOXP1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.156-10ng/mL
Sensitivity:0.062ng/mL
Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only