

## COL4A4 Rabbit pAb

Catalog No: #62063

Package Size: #62063-1 50ul #62063-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	COL4A4 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC IF ELISA
Species Reactivity	Human Mouse
Specificity	COL4A4 Polyclonal Antibody detects endogenous levels of COL4A4 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human Collagen IV alpha4.
Target Name	COL4A4
Other Names	COL4A4 ; Collagen alpha-4 ; IV ; chain
Uniprot	P53420
GeneID	1286
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

IHC 1:100-1:300; IF 1:200-1:1000; ELISA 1:5000

## Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using Collagen IV alpha4 Antibody. The picture on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of COS7 cells, using Collagen IV alpha4 Antibody. The picture on the right is blocked with the synthesized peptide.

---

Note: This product is for in vitro research use only