

Mouse Insulin Like Growth Factor 2 mRNA Binding Protein 3 (IGF2BP3) ELISA Kit



Catalog No: #EK20835

Orders: order@signalwayantibody.com

Package Size: #EK20835-1 48T #EK20835-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Insulin Like Growth Factor 2 mRNA Binding Protein 3 (IGF2BP3) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:6.35pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only