

Human Cadherin-related family member 2 (CDHR2) ELISA Kit

Catalog No: #EK20852

Package Size: #EK20852-1 48T #EK20852-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Cadherin-related family member 2 (CDHR2) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.062ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only