

# Human Large ribosomal subunit protein bL27m(MRPL27) ELISA Kit



Catalog No: #EK20845

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20845-1 48T #EK20845-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human Large ribosomal subunit protein bL27m(MRPL27) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range:0.781-50ng/mL

Sensitivity:0.35ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only