

Human Ran-binding protein 9(RANBP9) ELISA Kit

Catalog No: #EK20851



Package Size: #EK20851-1 48T #EK20851-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Ran-binding protein 9(RANBP9) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.123ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only