

HERC2 Rabbit pAb

Catalog No: #62066

Package Size: #62066-1 50ul #62066-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HERC2 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human, Mouse
Specificity	This antibody detects endogenous levels of human HERC2
Immunogen Description	Synthesized peptide derived from human HERC2 AA range: 646-696
Target Name	HERC2
Other Names	E3 ubiquitin-protein ligase HERC2 (HECT domain and RCC1-like domain-containing protein 2)
Uniprot	O95714;
GeneID	8924;
Calculated MW	532kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-15°C to -25°C/1 year(Do not lower than -25°C)

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only