

# Mouse Tumor Necrosis Factor Ligand Superfamily, Member 7 (TNFSF7) ELISA Kit



Catalog No: #EK20840

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20840-1 48T #EK20840-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Tumor Necrosis Factor Ligand Superfamily, Member 7 (TNFSF7) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range:7.81-500pg/mL

Sensitivity:3.56pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only