

Mouse Tumor Necrosis Factor Ligand Superfamily, Member 7 (TNFSF7) ELISA Kit



Catalog No: #EK20840

Orders: order@signalwayantibody.com

Package Size: #EK20840-1 48T #EK20840-2 96T

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Mouse Tumor Necrosis Factor Ligand Superfamily, Member 7 (TNFSF7) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Mouse |

Application Details

| |
|--|
| Detect Range:7.81-500pg/mL |
| Sensitivity:3.56pg/mL |
| Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids |
| Sample Volume: 50 μL |
| Assay Time:1-4.5h |
| Detection wavelength:450 nm |

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only