

Mouse Myelin Oligodendrocyte Glycoprotein (MOG) ELISA Kit



Catalog No: #EK20850

Orders: order@signalwayantibody.com

Package Size: #EK20850-1 48T #EK20850-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Myelin Oligodendrocyte Glycoprotein (MOG) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:10.7pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only