

SAK Antibody

Catalog No: #62068

Package Size: #62068-1 50ul #62068-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SAK Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Human; Mouse; Rat
Immunogen Description	Recombinant full length protein of human SAK
Target Name	PLK4
Uniprot	O00444
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

WB 1:500-1:2000

Images



Western blot analysis of SAK expression in MCF7 (A), HeLa (B), mouse spleen (C), rat brain (D) whole cell lysates. (Predicted band size: 104; 105; 108 kDa; Observed band size: 115 kDa)

Note: This product is for in vitro research use only