

Mouse Muscle-specific RING finger protein 1(Murf1) ELISA Kit



Catalog No: #EK20854

Orders: order@signalwayantibody.com

Package Size: #EK20854-1 48T #EK20854-2 96T

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mouse Muscle-specific RING finger protein 1(Murf1) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Mouse |

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.054ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

References

Detection Method:sandwich

Note: This product is for in vitro research use only