

Mouse Lactate Dehydrogenase D (LDHD) ELISA Kit

Catalog No: #EK20860



Package Size: #EK20860-1 48T #EK20860-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Lactate Dehydrogenase D (LDHD) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:3.12-200u/L

Sensitivity:1.3U/L

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only