

## CTHRC1 Rabbit pAb

Catalog No: #62082



Package Size: #62082-1 50ul #62082-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CTHRC1 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of CTHRC1 at Human, Mouse,Rat
Immunogen Description	Synthesized peptide derived from human CTHRC1
Target Name	CTHRC1 UNQ762/PRO1550
Uniprot	Q96CG8
GeneID	115908
Calculated MW	27kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:500-2000

Note: This product is for in vitro research use only