

# Rat Creatine Kinase,Brain (CKB) ELISA Kit

Catalog No: #EK20861



Package Size: #EK20861-1 48T #EK20861-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Creatine Kinase,Brain (CKB) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat

## Application Details

Detect Range:1.56-100ng/mL

Sensitivity:0.63ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only