

Human Arachidonate 12-lipoxygenase, 12R-type(ALOX12B) ELISA Kit



Catalog No: #EK20882

Orders: order@signalwayantibody.com

Package Size: #EK20882-1 48T #EK20882-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human Arachidonate 12-lipoxygenase, 12R-type(ALOX12B) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range	0.156-10ng/mL
Sensitivity	0.067ng/mL
Sample Type	Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids
Sample Volume	50 μL
Assay Time	1-4.5h
Detection wavelength	450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only