

# Human Lysine-specific demethylase 6B(KDM6B) ELISA Kit



Catalog No: #EK20870

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20870-1 48T #EK20870-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Human Lysine-specific demethylase 6B(KDM6B) ELISA Kit |
| Brief Description  | ELISA Kit   |
| Applications       | ELISA   |
| Species Reactivity | Human   |

## Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.112ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only