

Caspase-1 Mouse mAb

Catalog No: #62079

Package Size: #62079-1 50ul #62079-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | Caspase-1 Mouse mAb |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG2b,Kappa |
| Purification | Protein G |
| Applications | WB;IF;ELISA |
| Species Reactivity | Human;Mouse;Rat |
| Specificity | This antibody detects endogenous levels of Caspase-1 protein. |
| Immunogen Description | Synthesized peptide derived from human Caspase-1 |
| Target Name | CASP1 IL1BC IL1BCE |
| Uniprot | P29466 |
| GeneID | 834 |
| Calculated MW | 45kDa |
| Concentration | 1 mg/ml |
| Formulation | PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA |
| Storage | Store at -20°C/1 year |

Application Details

WB 1:500-2000; IF 1:100-500; ELISA 1:1000-5000

Images



Whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with anti-Caspase-1 antibody.

Note: This product is for in vitro research use only